

# PTPN14 antibody

Catalog Number: 114316

## Product name

PTPN14 antibody

WB: 1:200-1:2000

IHC: 1:20-1:200


## Specificity

Human, Mouse, Rat; other species not tested.

## Validations

## Antibody description

PTPN14 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, HeLa cells. Positive IHC detected in human placenta tissue, human brain tissue, human kidney tissue, human lung tissue, human ovary tissue, human skin tissue, human spleen tissue. Observed molecular weight by Western-blot: 150 kDa




## Preparation

This antibody was obtained by immunization of PTPN14 recombinant protein (Accession Number: XM\_047426370). Purification method: Antigen affinity purified.

K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:114316(PTPN14 antibody) at dilution of 1:400

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IHC

## Dilutions

Recommended Dilution:

Immunohistochemical of paraffin-embedded human placenta using Catalog No:114316(PTPN14 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded

# PTPN14 antibody



Catalog Number: 114316

---

human placenta using Catalog

No:114316(PTPN14 antibody) at dilution of 1:50  
(under 40x lens)