

TRIM41 antibody



Catalog Number: 116316

Product name

TRIM41 antibody

WB: 1:500-1:5000

IP: 1:200-1:2000

Specificity

Human, Mouse, Rat; other species not tested.


IHC: 1:20-1:200

IF: 1:10-1:100

Antibody description

TRIM41 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IHC detected in human skeletal muscle tissue. Positive WB detected in BxPC-3 cells, human placenta tissue, mouse liver tissue. Positive IP detected in human plasma tissue. Observed molecular weight by Western-blot: 70-72 kDa

Validations



Preparation

This antibody was obtained by immunization of TRIM41 recombinant protein (Accession Number: BC018765). Purification method: Antigen affinity purified.

BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:116316(TRIM41 antibody) at dilution of 1:800

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

ELISA, WB, IF, IHC, IP

Dilutions



Recommended Dilution:



Immunohistochemical of paraffin-embedded human skeletal muscle using Catalog No:116316(TRIM41 antibody) at dilution of 1:50 (under 10x lens)

TRIM41 antibody

Catalog Number: 116316



IgG (red).

Immunofluorescent analysis of HepG2 cells, using TRIM41 antibody Catalog No:116316 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit

IP Result of anti-TRIM41 (IP:Catalog No:116316, 4ug; Detection:Catalog No:116316 1:600) with human plasma tissue lysate 4000ug.