

Human MME/EPN (His tag) recombinant protein



Catalog Number: 518316

General Information

Gene Name Synonym

Nepriylsin; Atriopeptidase; Common acute lymphocytic leukemia antigen; CALLA; Enkephalinase; Neutral endopeptidase 24.11; NEP; Neutral endopeptidase; Skin fibroblast elastase; SFE; CD10; MME; EPN

Protein Construction

A DNA sequence encoding the Human MME/EPN Tyr52-Trp750 is expressed with a 8His tag at the N-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of PBS, pH7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human MME/EPN (His tag) recombinant protein