

Human LAG3 / CD223 / Lymphocyte activation gene 3 (Fc Tag) recombinant protein



Catalog Number: 501716

General Information

Gene Name Synonym

Protein FDC

Protein Construction

A DNA sequence encoding the human LAG3 (NP_002277.4) (Met1-Leu450) was expressed with the Fc region of human IgG1 at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> (3.0+91.0+3.1) % as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Val 29

Molecular Mass

The recombinant human LAG3 consists 660 amino acids and predicts a molecular mass of 72.3 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C.

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human LAG3 / CD223 / Lymphocyte activation gene 3 Protein (Fc Tag)