

# Human CD3d/CD3 delta (Fc & FLAG Tag) recombinant protein



Catalog Number: 504716

## General Information

### Protein Construction

A DNA sequence encoding the human CD3D (Met1-Ala105) was fused with the C-terminal flag-tagged Fc region of human IgG1 at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

> 85 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Phe 22

### Molecular Mass

The recombinant human CD3D consists of 335 amino acids and predicts a molecular mass of 37.8

KDa. It migrates as an approximately 47-51 KDa band in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human CD3d / CD3 delta Protein (Fc & FLAG Tag)  
SDS-PAGE