

STK31 antibody

Catalog Number: 115716

Product name

STK31 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

STK31 Rabbit Polyclonal antibody. Positive IP detected in mouse testis tissue. Positive WB detected in mouse testis tissue, rat testis tissue. Observed molecular weight by Western-blot: 115 kDa

Preparation

This antibody was obtained by immunization of STK31 recombinant protein (Accession Number: NM_031414). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:500-1:5000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:115716(STK31 Antibody) at dilution of 1:600



IP Result of anti-STK31 (IP:Catalog No:115716, 4ug; Detection:Catalog No:115716 1:1000) with mouse testis tissue lysate 2400ug.