

Anti-RCA antibody

Catalog Number: 176716

Product name

Anti-RCA antibody

Specificity

Arabidopsis thaliana, Nicotiana tabacum, Oryza sativa, Zea mays, Spinacia oleracea, Morus alba

Antibody description

Rabbit polyclonal antibody to RCA

Preparation

This antigen of this antibody was synthetic peptide within c-terminal arabidopsis thaliana rubisco activase.

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 50% Glycerol. Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:10,000-1:50,000

IHC-P: 1:50-1:200

Validations




Fig1: Western blot analysis on *A. thaliana* tissue lysates using anti-Rubisco activase rabbit polyclonal antibody at 1/10,000 dilution.




Fig2: Immunohistochemical analysis of paraffin-embedded *A. thaliana* using anti-Rubisco activase rabbit polyclonal antibody.




Fig3: Immunohistochemical analysis of paraffin-embedded *A. thaliana* using anti-Rubisco activase rabbit polyclonal antibody.