

# Cystatin F antibody



Catalog Number: 109716

## Product name

Cystatin F antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

Cystatin F Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells, Jurkat cells, mouse uterus tissue. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 19 kDa

## Preparation

This antibody was obtained by immunization of Cystatin F recombinant protein (Accession Number: BC015507). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:109716(CST7 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:109716(CST7 antibody) at dilution of 1:50 (under 10x lens)