

Cystatin F antibody



Catalog Number: 109716

Product name

Cystatin F antibody

ELISA, WB, IHC

Specificity

Human, Mouse; other species not tested.

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

IHC: 1:20-1:200

Antibody description

Cystatin F Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells, Jurkat cells, mouse uterus tissue. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 19 kDa

Validations

Preparation

This antibody was obtained by immunization of Cystatin F recombinant protein (Accession Number: BC015507). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Polyclonal


Ig Type

Rabbit IgG

Applications



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:109716(CST7 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:109716(CST7 antibody) at dilution of 1:50 (under 10x lens)