

# NCALD antibody



Catalog Number: 113026

## Product name

NCALD antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

NCALD Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, human brain tissue, human cerebellum tissue, human heart tissue, mouse cerebellum tissue, rat brain tissue. Observed molecular weight by Western-blot: 18-22 kDa

## Preparation

This antibody was obtained by immunization of NCALD recombinant protein (Accession Number: XM\_017013902). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions


Recommended Dilution:

WB: 1:1000-1:10000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113026(NCALD antibody) at dilution of 1:2000



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113026(NCALD Antibody) at dilution of 1:2000