

Spesolimab reference antibody



Catalog Number: 177026

Product name

Spesolimab reference antibody

Specificity

Human IL-1RL2 / IL-36R

Antibody description

Anti-IL-1RL2 / IL-36R Antibody (Spesolimab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

1×PBS, pH 6.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG

Applications


ELISA

Validations



SDS-PAGE

Anti-IL-1RL2 / IL-36R Antibody (spesolimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.



SEC-HPLC

The purity of Anti-IL-1RL2 / IL-36R Antibody (spesolimab) is more than 95%, determined by SEC-HPLC.