

Human ERP44/KIAA0573/TXNDC4 (His tag) recombinant protein



Catalog Number: 519026

General Information

Gene Name Synonym

Thioredoxin domain-containing protein 4, ER protein 44, KIAA0573, TXNDC4

Protein Construction

A DNA sequence encoding the Human ERP44/KIAA0573/TXNDC4 Glu30-Asp402 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM Tris HCl \square 10% glycerol, PH7.5.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human ERP44/KIAA0573/TXNDC4 (His tag)
recombinant protein