

# HLA-DRB1 antibody



Catalog Number: 111426

## Product name

HLA-DRB1 antibody

## Specificity

Human; other species not tested.

## Antibody description

HLA-DRB1 Rabbit Polyclonal antibody. Positive WB detected in human spleen tissue, A375 cells, Raji cells. Positive IP detected in Raji cells. Observed molecular weight by Western-blot: 27-30 kDa

## Preparation

This antibody was obtained by immunization of HLA-DRB1 recombinant protein (Accession Number: XM\_047444769). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



human spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:111426(HLA-DRB1 Antibody) at dilution of



IP Result of anti-HLA-DRB1 (IP:Catalog No:111426, 4ug; Detection:Catalog No:111426 1:500) with Raji cells lysate 2000ug.