

# Human CD16a/FCGR3A (176 Phe, His Tag) recombinant protein



Catalog Number: 503426

## General Information

### Gene Name Synonym

CD16a antigen; Fc-gamma RIII-alpha; FcR-10; IgG Fc receptor III-2

### Protein Construction

A DNA sequence encoding the extracellular domain (Met 1-Gln 208) of human CD16a (P08637-1) precursor was expressed with a C-terminal polyhistidine tag. It is identical to FCGR3A158F/V in the reference.

### Organism

Human

### Expression Host

CHO Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Gly 17

## Molecular Mass

The secreted recombinant human CD16a consists of 203 amino acids and predicts a molecular mass of 23.3 KDa.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human CD16a / FCGR3A Protein (His Tag) SDS-PAGE