

Anti-XCL1 antibody

Catalog Number: 176426

Product name

Anti-XCL1 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to XCL1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human xcl1/xcl2. aa range:51-100

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.; Western Blot analysis of 293 cells using Lymphotactin Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at




Fig2.; Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:100