

# Anti-Lysozyme Rabbit antibody



Catalog Number: 167426

## Product name

Anti-Lysozyme Rabbit antibody

WB: 1/1000

ICC/IF: 1/50

## Specificity

Human

IP: 1/20

## Antibody description

Lysozyme Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human Lysozyme

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications


WB, ICC/IF, IP

## Dilutions

## Validations



Western blot detection of Lysozyme in HL-60 cell lysates using Lysozyme Rabbit pAb(1:1000 diluted). Predicted band size:17KDa. Observed



band size:17KDa.

Western blot detection of Lysozyme in Jurkat,C6,3T3 cell lysates using Lysozyme Rabbit pAb(1:1000 diluted). Predicted band size:17KDa. Observed band size:17KDa.