

MLANA antibody

Catalog Number: 107426

Product name

MLANA antibody

Specificity

Human; other species not tested.

Antibody description

MLANA Mouse Monoclonal antibody. Positive IHC detected in human malignant melanoma tissue.

Preparation

This antibody was obtained by immunization of MLANA recombinant protein (Accession Number: NM_005511). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:50-1:500

Validations



Immunohistochemistry of paraffin-embedded human malignant melanoma tissue slide using Catalog No:107426(MLANA Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemistry of paraffin-embedded human malignant melanoma tissue slide using Catalog No:107426(MLANA Antibody) at dilution of 1:200 (under 40x lens).