

Human FAM3D / Oit1 (His Tag) recombinant protein



Catalog Number: 502826

General Information

Protein Construction

A DNA sequence encoding the human FAM3D (Q96BQ1-1) (Met 1-Phe 224) was fused with a polyhistidine tag at the C-terminus

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Tyr 26

Molecular Mass

The recombinant human FAM3D consists of 210 amino acids and has a predicted molecular mass of 23.4 kDa. In SDS-PAGE under reducing

conditions, rhFAM3D migrates as an approximately 27 kDa band due to glycosylation

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human FAM3D Protein (His Tag) SDS-PAGE