

Anti-Eif4ebp1 antibody



Catalog Number: 175826

Product name

Polyclonal

Anti-Eif4ebp1 antibody

Ig Type

Specificity

Rabbit IgG

Mouse, Rat, Monkey

Applications

Antibody description

WB, IHC-P

Rabbit polyclonal antibody to Eif4ebp1

Dilutions

Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from mouse 4e-bp1. aa range:35-84

IHC-P: 1:100-1:300

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Validations




Fig1.; Western Blot analysis of various cells using 4E-BP1 Polyclonal Antibody diluted at 1: 500