

Anti-Eif4ebp1 antibody



Catalog Number: 175826

Product name

Anti-Eif4ebp1 antibody

Specificity

Mouse, Rat, Monkey

Antibody description

Rabbit polyclonal antibody to Eif4ebp1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from mouse 4e-bp1. aa range:35-84

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.: Western Blot analysis of various cells using 4E-BP1 Polyclonal Antibody diluted at 1: 500