

Amivantamab reference antibody



Catalog Number: 176826

Product name

Amivantamab reference antibody

Specificity

Human HGFR / c-Met

Antibody description

Anti-HGFR / c-Met Antibody (Amivantamab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG1

Applications


ELISA

Validations




SDS-PAGE

Anti-HGFR / c-Met Antibody (amivantamab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



SEC-HPLC

The purity of Anti-HGFR / c-Met Antibody (amivantamab) is more than 95%, determined by



SEC-HPLC.

Bioactivity: ELISA

Immobilized human cMet His at 2 ug/mL can bind Anti-HGFR / c-Met Antibody (amivantamab), EC50=0.006091 ug/mL.