

# Anti-PCTAIRE1 Rabbit antibody



Catalog Number: 167826

## Product name

Anti-PCTAIRE1 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

PCTAIRE1 Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human PCTAIRE1

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of PCTAIRE1 in A549,HL-60,U251 cell lysates using PCTAIRE1 Rabbit mAb(1:1000 diluted).Predicted band size:56kDa.Observed band size:65kDa.