

ARHGEF2 Polyclonal Antibody



Catalog Number: 161126

Product name

ARHGEF2 Polyclonal Antibody

Specificity

Mouse

Antibody description

Polyclonal antibody to ARHGEF2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-350 of human ARHGEF2 (NP_004714.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions

IHC 1:50 - 1:100

Validations



Immunohistochemistry - ARHGEF2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse lung using ARHGEF2 antibody at dilution of 1:100 (40x lens).