

# SERCA2 antibody

Catalog Number: 115126

---

## Product name

SERCA2 antibody

## Immunogen

[Human SERCA2,ATP2A2 Recombinant protein \(His tag\)](#)

## Specificity

Human,Mouse,Rat; other species not tested.

## Antibody description

SERCA2 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue. Observed molecular weight by Western-blot: 100kd

## Preparation

This antibody was obtained by immunization of SERCA2 recombinant protein (Accession Number: NM\_170665). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:1000-1:10000

## Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:115126(SERCA2,ATP2A2 antibody) at dilution of 1:3000