

# Anti-JAK1 (7B1) Mouse antibody



Catalog Number: 168126

## Product name

WB, IP

Anti-JAK1 (7B1) Mouse antibody

## Dilutions

## Specificity


WB: 1/500

Human, Rat

## Antibody description

## Validations

JAK1 (7B1) Mouse monoclonal antibody



## Preparation

Antigen: Purified recombinant human Jak1 protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Immunoprecipitation analysis of HeLa cell lysates using Jak1 mouse mAb (168126). Jak1 mouse mAb (200622-8B8) was used for the western blot

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications



Western blot analysis of extracts from PC12,C6,Ramos and Jurkat cell lysates using Jak1 mouse mAb (1:500 diluted).Predicted band size:130kDa.Observed band size:130kDa.