

Anti-INS antibody



Catalog Number: 175526

Product name

Anti-INS antibody

Specificity

Human, Mouse, Rat, Pig

Antibody description

Rabbit polyclonal antibody to INS

Preparation

This antigen of this antibody was porcine pancreas \square full length

Formulation

Liquid, 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

Storage

Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB:1:500-2000

Validations




Fig1: Sample;; Lane 1:Islet A lysates at 20ug;; line2,Islet B lysates at 20ug;; Primary: Anti-Insulin at 1:300 dilution;; Secondary: HRP conjugated Goat-Anti-rabbit IgG(bs-0295G-HRP) at 1: 5000 dilution;; Predicted band size:5.8 kD Observed band size:5.8 kD