

UBE2G2 antibody

Catalog Number: 116526

Product name

UBE2G2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

UBE2G2 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Positive IHC detected in human prostate cancer tissue.

Preparation

This antibody was obtained by immunization of UBE2G2 recombinant protein (Accession Number: NM_003343). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:116526(UBE2G2 antibody) at dilution of 1:100 (under 25x lens)



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116526(UBE2G2 Antibody) at dilution of 1:1000