

# KCNB1-Specific antibody



Catalog Number: 111926

---

## Product name

KCNB1-Specific antibody

## Specificity

Human; other species not tested.

## Antibody description

KCNB1-Specific Rabbit Polyclonal antibody.  
Positive IF detected in SH-SY5Y cells, HEK-293 cells.

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_004975).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IF

## Dilutions

Recommended Dilution:

IF: 1:10-1:100

## Validations



Immunofluorescent analysis of SH-SY5Y cells using Catalog No:111926(KCNB1-Specific Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG