

MYF6 antibody



Catalog Number: 112926

Product name

MYF6 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MYF6 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 27 kDa

Preparation

This antibody was obtained by immunization of MYF6 recombinant protein (Accession Number: NM_002469). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:112926(MYF6 antibody) at dilution of 1:500