

RRN3 antibody



Catalog Number: 114926

Product name

RRN3 antibody

Immunogen

[Human RRN3 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

RRN3 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, Jurkat cells, K-562 cells, mouse testis tissue. Observed molecular weight by Western-blot: 74 kDa

Preparation

This antibody was obtained by immunization of RRN3 recombinant protein (Accession Number: NM_018427). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114926(RRN3 Antibody) at dilution of 1:600