

# Golimumab reference antibody



Catalog Number: 176926

## Product name

Golimumab reference antibody

## Specificity

Human TNFSF2 / TNF $\alpha$

## Antibody description

Anti-TNFSF2 / TNF $\alpha$  Antibody (Golimumab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

1 $\times$ PBS, pH 6.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG

## Applications


## ELISA

## Validations



## SDS-PAGE

Anti-TNFSF2 / TNF $\alpha$  Antibody (golimumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.



## SEC-HPLC

The purity of Anti-TNFSF2 / TNF $\alpha$  Antibody (golimumab) is more than 95%, determined by SEC-HPLC.