

AK5 antibody

Catalog Number: 107926

Product name

AK5 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

AK5 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 63kd

Preparation

This antibody was obtained by immunization of AK5 recombinant protein (Accession Number: BC012467). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107926(AK5 antibody) at dilution of 1:300