

EHD4 antibody



Catalog Number: 110226

Product name

EHD4 antibody

Immunogen

[Human EHD4 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

EHD4 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in mouse testis tissue, HeLa cells. Observed molecular weight by Western-blot: 61 kDa

Preparation

This antibody was obtained by immunization of EHD4 recombinant protein (Accession Number: NM_139265). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:110226(EHD4 antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using EHD4 antibody Catalog No:110226 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).