

# LILRB5 antibody



Catalog Number: 112226

---

## Product name

LILRB5 antibody

## Specificity

human; other species not tested.

## Antibody description

LILRB5 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 64kd

## Preparation

This antibody was obtained by immunization of LILRB5 recombinant protein (Accession Number: NM\_006840). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:112226(LILRB5 antibody) at dilution of 1:1500