

SIGLEC7 antibody

Catalog Number: 115226

Product name

SIGLEC7 antibody

WB: 1:500-1:5000

IHC: 1:20-1:200


Specificity

Human; other species not tested.

Validations

Antibody description

SIGLEC7 Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue, human liver tissue, human placenta tissue. Positive WB detected in U-937 cells, A549 cells, human liver tissue, human placenta tissue, K-562 cells, L02 cells. Observed molecular weight by Western-blot: 75 kDa




Preparation

This antibody was obtained by immunization of SIGLEC7 recombinant protein (Accession Number: NM_014385). Purification method: Antigen affinity purified.

U-937 cells were subjected to SDS PAGE followed by western blot with Catalog No:115226(SIGLEC7 antibody) at dilution of 1:1000

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:

Immunohistochemistry of paraffin-embedded human tonsillitis slide using Catalog No:115226(SIGLEC7 Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human tonsillitis slide using Catalog No:115226(SIGLEC7 Antibody) at dilution of 1:50

SIGLEC7 antibody



Catalog Number: 115226

No:115226(SIGLEC7 Antibody) at dilution of 1:50