

# TRA2A antibody

Catalog Number: 116226

---

## Product name

TRA2A antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TRA2A Rabbit Polyclonal antibody. Positive WB detected in A549 cells. Observed molecular weight by Western-blot: 33 kDa

## Preparation

This antibody was obtained by immunization of TRA2A recombinant protein (Accession Number: NM\_013293). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:116226(TRA2A antibody) at dilution of 1:3000