

# CYP2D6 antibody



Catalog Number: 107226

## Product name

CYP2D6 antibody

## Immunogen

[Human CYP2D6 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

CYP2D6 Mouse Monoclonal antibody. Positive WB detected in HepG2 cells, fetal human brain tissue. Observed molecular weight by Western-blot: 56 kDa

## Preparation

This antibody was obtained by immunization of CYP2D6 recombinant protein (Accession Number: NM\_000106). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:107226(CYP2D6 Antibody) at dilution of 1:1000