

# FFAR2 antibody

Catalog Number: 110626

---

## Product name

FFAR2 antibody

## Specificity

Human; other species not tested.

## Antibody description

FFAR2 Rabbit Polyclonal antibody. Positive IHC detected in human spleen tissue, human colon tissue.

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_005306).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:50-1:500

## Validations



Immunohistochemistry of paraffin-embedded human spleen tissue slide using Catalog No:110626(FFAR2 Antibody) at dilution of 1:100 (under 10x lens).



Immunohistochemistry of paraffin-embedded human spleen tissue slide using Catalog No:110626(FFAR2 Antibody) at dilution of 1:100 (under 40x lens).