

# MIP antibody



Catalog Number: 112626

## Product name

MIP antibody

## Immunogen

[Human MIP Recombinant protein \(GST tag\)](#)

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

MIP Rabbit Polyclonal antibody. Positive WB detected in mouse eye tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 26-28kd

## Preparation

This antibody was obtained by immunization of MIP recombinant protein (Accession Number: NM\_012064). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



mouse eye tissue were subjected to SDS PAGE followed by western blot with Catalog No:112626(MIP antibody) at dilution of 1:500