

ST3GAL6 antibody



Catalog Number: 115626

Product name

ST3GAL6 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

ST3GAL6 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue, human heart tissue. Observed molecular weight by Western-blot: 50-55kd

Preparation

This antibody was obtained by immunization of ST3GAL6 recombinant protein (Accession Number: NM_001323368). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:115626(ST3GAL6 antibody) at dilution of 1:200