

LIF Polyclonal Antibody



Catalog Number: 161326

Product name

LIF Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to LIF

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 23-202 of human LIF (NP_002300.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions

IHC 1:50 - 1:100

Validations



Immunohistochemistry - LIF Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human prostate using LIF antibody at dilution of 1:100 (40x lens).