

# Human MRPL44/L44mt (His Tag) recombinant protein



Catalog Number: 501326

## General Information

### Protein Construction

A DNA sequence encoding the human MRPL44 (Q9H9J2) (Met1-Ser332) was expressed with a C-terminal polyhistidine tag.

### Organism

Human

### Expression Host

Baculovirus-Insect Cells

## QC Testing

### Purity

> 85 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Val 31

### Molecular Mass

The secreted recombinant human MRPL44 consists of 312 amino acids and predicts a molecular mass of 35.7 KDa. The apparent

molecular mass of the protein is approximately 40 KDa in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile 20mM, Tris 500mM NaCl, pH 8.0.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human MRPL44 Protein (His Tag) SDS-PAGE