

H3N2 HA (His Tag) recombinant protein



Catalog Number: 502326

General Information

Protein Construction

A DNA sequence encoding the Influenza A virus (A/Guangdong-Luohu/1256/2009(H3N2)) hemagglutinin (AFM72872.1) (Met1-Trp530), termed as HA, was expressed with a polyhistidine tag at the C-terminus.

Organism

H3N2

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 85 % as determined by SDS-PAGE.

Endotoxin

< 1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 17

Molecular Mass

The recombinant hemagglutinin of Influenza A virus (A/Guangdong-Luohu/1256/2009(H3N2)) consists 525 amino acids and predicts a molecular

mass of 59.2 kDa.

Formulation

- Lyophilized from sterile 20 mM Tris, 500 mM NaCl, 10 % glycerol, pH 7.0.
- 5 % trehalose and mannitol are added as protectants before lyophilization.
 - Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H3N2 (A/Guangdong-Luohu/1256/2009) Hemagglutinin / HA Protein (His Tag)