

LRR 1 antibody



Catalog Number: 112326

Product name

LRR 1 antibody

Immunogen

[Human PPIL5 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

LRR 1 Rabbit Polyclonal antibody. Positive WB detected in U-937 cells. Positive IHC detected in human pancreas cancer tissue. Observed molecular weight by Western-blot: 47kd

Preparation

This antibody was obtained by immunization of LRR 1 recombinant protein (Accession Number: NM_203467). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



U-937 cells were subjected to SDS PAGE followed by western blot with Catalog No:112326(PPIL5 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:112326(PPIL5 antibody) at dilution of 1:50 (under 10x lens)