

# Caspase 8 Polyclonal Antibody



Catalog Number: 160726

## Product name

Caspase 8 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to Caspase 8

## Preparation

Antigen: Recombinant protein of mouse Caspase 8

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

## Validations



Western blot - Caspase 8 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using Caspase 8 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - Caspase 8 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using Caspase 8 antibody at dilution of 1:100 (40x lens).