

Anti-DGKZ/DGK-zeta Rabbit antibody



Catalog Number: 167726

Product name

Anti-DGKZ/DGK-zeta Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

DGKZ/DGK-zeta Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human DGKZ/DGK-zeta

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, FC, IP

Dilutions


WB: 1/1000

IHC: 1/500

FC: 1/50

IP: 1/20

Validations



Western blot detection of DGKZ/DGK-zeta in HeLa,A549,HL-60,U2OS,U87-MG cell lysates using DGKZ/DGK-zeta Rabbit mAb(1:1000 diluted).Predicted band size:104kDa.Observed band size:124,114kDa.