

C3orf43 antibody



Catalog Number: 108726

Product name

C3orf43 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

C3orf43 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells, HeLa cells. Positive WB detected in mouse heart tissue, human plasma tissue. Observed molecular weight by Western-blot: 25-30 kDa

Preparation

This antibody was obtained by immunization of C3orf43 recombinant protein (Accession Number: NM_001077657). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:108726(C3orf43 Antibody) at dilution of 1:600



Immunofluorescent analysis of HepG2 cells using Catalog No:108726(C3orf43 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG