

Human PKM/OIP3/PK2/PK3/PKM2 (His tag) recombinant protein



Catalog Number: 518726

General Information

Gene Name Synonym

Pyruvate kinase PKM, CTHBP, OIP-3, THBP1

Protein Construction

A DNA sequence encoding the Human PKM/OIP3/PK2/PK3/PKM2 Ser2-Pro531 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of PBS, pH7.0, 10% glycerol.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human PKM/OIP3/PK2/PK3/PKM2 (His tag)
recombinant protein