

# CCBP2 antibody



Catalog Number: 107036

## Product name

CCBP2 antibody

## Specificity

Human; other species not tested.

## Antibody description

CCBP2 Mouse Monoclonal antibody. Positive IHC detected in human tonsillitis tissue.

## Preparation

This antibody was obtained by immunization of CCBP2 recombinant protein (Accession Number: NM\_001296). Purification method: Caprylic acid/ammonium sulfate precipitation.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:107036(CCBP2 Antibody) at dilution of 1:400 (under 10x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:107036(CCBP2 Antibody) at dilution of 1:400 (under 40x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).