

MAPK7 antibody



Catalog Number: 112436

Product name

MAPK7 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MAPK7 Rabbit Polyclonal antibody. Positive WB detected in SKOV-3 cells, HEK-293 cells, mouse heart tissue, mouse testis tissue. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of MAPK7 recombinant protein (Accession Number: NM_002749). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



SKOV-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:112436(MAPK7 Antibody) at dilution of 1:600